



MiniCute imager

Small, Light, and Easy

- **Switchable filters**

Switchable UV and blue light filters

- **Compact:**

Smaller dimension, save your working area

- **Convenient:**

WiFi data transfer, easy to download your data and print via WiFi

- **Easy-to-use:**

No complex interface, just operate digital camera

- **Open system:**

High compatibility to other branded transilluminators



MiniCute imager is an easy gel image system , recommend to use with UltraSlim® LED Illuminator, UltraBright LED Transilluminator and UltraBright UV Transilluminator



MiniCute imager, compact size, compatible with Maestrogen's UV and blue LED transilluminators.

- 100% fit to all Maestrogen's LED or UV illuminators
- Branded digital camera, 100% dark room for photography
- Light is convenient to move around (option)
- adaptor + fake battery for DSC (option)

Order Information and Specification:

Products	miniCute imager
Cat No.	MC-02
Outer dimensions (mm)	275 (W) x 238(D) x 250(H)
Filter	Amber filter and UV-cut filter
Power	5V/5.5W
Weight (Kg)	1.1kg
Digital camera:	(1) Effective pixel : 12.2 M (2) Lens: Fujinon 15x lens, F3.1 (Wide) - F5.6 (Telephoto) (3) Mode: Program (P), Aperture Priority (A), Shutter Priority (S), Manual (M) (4) ISO : Auto/100/200/400/800/1600/3200/6400 (5) Video Output: NTSC / PAL selectable (6) Storage Media: SD/SDHC memory card/ WiFi Transfer (7) Storage Format: JPEG (ExifVer 2.2)

Related product:

Products	UltraSlim® LED Illuminator	UltraBright LED Transilluminator	UltraBright UV Transilluminator	adaptor + fake battery for DSC (option)
Cat No.	SLB-01/SLB-01W	LB-16	MLB-21/MB-21 MLB-16/MB-16	MCD-01
Description	viewing stained gels on the laboratory bench, or within gel documentation systems			Power supply for DSC