



Now including a full range of
microcentrifuge tubes, caps and closures










Fisherbrand[®] SureOne Pipet Tips



Fisherbrand[™] SureOne[™] pipet tips are a comprehensive line of universal fit pipet tips offered in many packaging configurations. Tips are available in volume range from 5 µl to 10 ml.

- Sterile product certified free of RNase, DNase, Pyrogen, Bioburden, PCR Inhibitors, Endotoxin and DNA
- Genomic tips and environmentally-friendly reload packaging also available
- Non sterile product certified free of RNase, DNase, DNA and Pyrogen
- Optimal fit with the Fisherbrand pipettes

SureOne Filter Tips

Cat. No.	Tip Volume	Packaging	Color	
➤ 11917734	0.1 - 5 µl	Filtered Sterile Filter Sterile Reload	Clear	 Universal Fits, Micropoint Nano tip, graduated at 2.5 µl. Tip length 1.22 in, 3.10 cm
➤ 11903466	0.1 - 10 µl	Filtered Sterile Filter Sterile Reload	Clear	 Universal Fits, Micropoint tip, graduated at 2 µl. Tip length 1.23 in, 3.13 cm
➤ 11913466	0.1 - 10 µl	Filtered Sterile Filter Sterile Reload	Clear	 Universal Fits, Micropoint tip graduated at 2, 5, 10, 20 µl, Extended length. Tip length 1.80 in, 4.57 cm
➤ 11977714	0.1 - 20 µl	Filtered Sterile *Filter Sterile Reload	Clear	 Universal Fits, Micropoint tip graduated at 2, 5, 10, 20 µl, Extended length. Tip length 1.80 in, 4.57 cm
➤ 11943466	2 - 20 µl	Filtered Sterile **Filter Sterile Reload	Clear	 Universal Fits, Beveled tip. Tip length 2.01 in, 5.11 cm
➤ 11953466	10 - 100 µl	Filtered Sterile Filter Sterile Reload	Clear	 Universal Fits, Beveled tip. Tip length 2.01 in, 5.11 cm
➤ 11963466	20 - 200 µl	Filtered Sterile Filter Sterile Reload	Clear	 Universal Fits, Micropoint tip graduated at 10, 50, 100, 200 µl. Tip length 2.30 in, 5.84 cm
➤ 11987714	30 - 300 µl	Filtered Sterile Filter Sterile Reload	Clear	 Universal Fits, Micropoint tip graduated at 10, 50, 100, 200 µl. Tip length 2.30 in, 5.84 cm
➤ 11973466	100 - 1000 µl	Filtered Sterile Filter Sterile Reload	Clear	 Universal Fits, Micropoint tip graduated at 100, 200, 500, 1000 µl. Tip length 3.50 in, 8.89 cm
➤ 11927754	1 - 5 ml	Filtered Sterile, Bulk	Clear	 Fits Gilson Pipetman and Gilson Style pipettes. Tip length 4.60 in, 11.70 cm

*Use empty rack 10 - 20 µl: 02-707-445

**Use empty rack 20 - 300 µl: 02-707-444

➤ **New!**

DISTRIBUIDOR AUTORIZADO
LABOTAQ, S.L.
TELF 954 31 72 16
PEDIDOS@LABOTAQ.COM
WWW.LABOTAQ.COM

 **labotaq**
molecular biology specialist